## **Biochemichal Conversion**

The department of Biochemical Conversion focuses on research for the provision of energy from biomass with the help of microorganisms. Here, the main focus lies on the technology used for biogas production and utilization.



12 L lab fermenters

To ensure the quality of the manifold research topics, the department is divided into the following work groups:

- Characterization and development of anaerobic processes
- System optimization •
- Process monitoring and simulation •
- Biogas technology

Additionally, the work group "Microbiology of Anaerobic Systems" (MicAS) works in the Biochemical Conversion department - in close collaboration with the Helmholtz Centre for Environmental Research (UFZ).

With support from



by decision of the German Bundestag

### Contact

Head of Department "Biochemical Conversion" Dr. agr. Peter Kornatz E-Mail: peter.kornatz@dbfz.de

#### **Biogas lab:**

Hed of working group "CAP" Dr. Nils Engler E-Mail: nils.engler@dbfz.de

Technical management Biogas lab Katrin Strach E-Mail: katrin.strach@dbfz.de

**Research biogas plant:** Head of working group "Biogas technology" Dr. Marcel Pohl E-Mail: marcel.pohl@dbfz.de

Technical management Research biogas plant Christian Krebs E-Mail: christian.krebs@dbfz.de

## **Further Information**

For additional information about the "Deutsches Biomasseforschungszentrum" and ongoing as well as completed projects, please refer to www.dbfz.de.

Further information about the Biochemical Conversion department and the biogas research plant can be obtained by scanning the QR code below.



Deutsches Biomasseforschungszentrum DBFZ gemeinnützige GmbH



**Applied Research in the Department** of Biochemical Conversion



Technical facilities of biogas lab and biogas research plant

### Biogas research plant

### Facilities

- 2 main fermenters (190 m<sup>3</sup>)
- 1 plug-flow fermenter (53 m<sup>3</sup>)
- 1 small-sized fermenter (88 m<sup>3</sup>)
- 1 secondary fermenter (215 m<sup>3</sup>)
- 1 manure storage tank (174 m<sup>3</sup>)

**Biogas** laboratory

Methane forming potential

Determination of alcohols

Trace gases (CH<sub>4</sub>, N<sub>2</sub>O, NH<sub>2</sub>)

Emissions from biogas plants

Sugar analytics

**Emission measurements:** 

۲

٠

۲

•

۲

•

•

Analytics of biogas feedstocks and digestates:

Composition analysis (Weender, van Soest)

Microscopy and physiological characterization

Emissions from digestate (field measurements)

- 1 digestate storage tank (215 m<sup>3</sup>)
- CHP (pilot injection gas engine) with electric power of 75 kW output, equipped with a flare



Standard analytics (TS, VS, TKN, NH<sub>4</sub>N, Titration of VOA, pH, etc.)

### Areas of research

- Process stabilization and mass balances and efficiency increase
- Scale-up from lab to industrial scale
- Use of microbial and abiotic parameters for process monitoring and control
- Flexibilization of biogas production
- Investigation of mixing techniques & mixing properties
- New feedstocks for anaerobic digestion
- Extreme process conditions
- New applications of sensors and moduls for process monitoring and control



## Areas of research

- Improvement of the efficiency of the biogas process
- Process optimization and process development for special substrates
- Effects of additive addition (trace elements, nutrients and enzymes)
- Comparison of different technologies for biogas production
- Effects of feedstock pretreatment
- Microbiology (in cooperation with UFZ)
- Development of analysis methods
- Emission measurements

## Technical capabilities

- Mesophilic and thermophilic operation
- Process control via a process control system
- Various programmable modes of operation
- Online measurements (e.g. pH, temperature)
- Determination of biogas quality and quantity for each fermenter
- Various options for sampling and technical modification
- Parallel operations due to the identical construction of the main fermenters
- Fittings for additional gas utilization technologies or biogas treatment



# Technical facilities

#### Batch assays:

• 70 eudiometers, 60 AMPTS

#### **Continuous fermentation tests:**

- 46 CSTR up to 15 L
- 2 CSTR up to 35 L
- 6 CSTR up to 45 L
- 1 CSTR up to 500 L
- 1 IC reactor 23 L
- 1 plug-flow fermenter 180 L

